

**TUBERCULOCIDAL ACTIVITY:** This product exhibits disinfectant efficacy against *Mycobacterium bovis* (BCG) in 3 minutes at 20° C when used as directed on hard, nonporous, inanimate surfaces in the presence of 5% blood serum.

**BACTERICIDAL AND \*VIRUCIDAL ACTIVITY:** This product exhibits disinfectant efficacy in 3 minutes against the following organisms when used as directed on hard, non-porous, inanimate surfaces in the presence of 5% blood serum: *Pseudomonas aeruginosa* (ATCC 15442), *Staphylococcus aureus* (ATCC 6538), *Salmonella enterica* (ATCC 10708), Methicillin-Resistant *Staphylococcus aureus* (MRSA, ATCC 33592), Influenza A2 Virus (A/2/Japan/305/57, strain H2N2), Human Immunodeficiency Virus Type 1 (HIV-1) (strain IIIB), Herpes Simplex Virus Type 2 (strain G, ATCC VR-734), Human Coronavirus (ATCC VR-740, Strain 229E), Feline Calicivirus (ATCC VR-782, Strain F-9), Norovirus (Norwalk virus), and Rhinovirus type 37 (ATCC VR-1147, Strain 151-1).  
**\*KILLS HIV ON PRE-CLEANED HARD, NON-POROUS SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health care settings or other settings in which there is an expected likelihood of soiling of hard, non-porous surfaces/objects with blood or body fluids and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1.

**SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 ON SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS.** **Personal Protection:** When handling items soiled with blood or body fluids, use disposable latex gloves, gowns, masks and eye protection. **Cleaning Procedure:** Blood and other body fluids must be thoroughly cleaned from hard, non-porous surfaces and objects before application of this product. **Disposal of Infectious Materials:** Blood and other body fluids, cleaning materials and clothing must be autoclaved and disposed of according to federal, state and local regulations for infectious waste disposal. **Contact Time:** Allow surfaces to remain wet for 3 minutes.

Disinfects and controls odors on hard, nonporous surfaces such as metal surfaces, stainless steel surfaces, glazed porcelain, glazed ceramic tile, plastic surfaces, bathroom fixtures, shower stalls, empty garbage cans and diaper pails, furniture in sick rooms and nurseries, fishing equipment and sports equipment. USED IN: HOSPITALS AND OTHER MEDICAL, INDUSTRIAL, INSTITUTIONAL, COMMERCIAL AND RESIDENTIAL FACILITIES. **DIRECTIONS FOR USE:** It is a violation of federal law to use this product in a manner inconsistent with its labeling. This product is not for use on critical or semi-critical medical device surfaces.

HOLD CONTAINER UPRIGHT. SHAKE OCCASIONALLY DURING USE. **FOR TUBERCULOCIDAL USE:** Remove visible soil and spray from 6-8" until thoroughly wet. Allow surface to remain undisturbed for at least 3 minutes at 20° C before wiping.

**FOR DISINFECTING and deodorizing:** remove visible soil and spray from 6-8" until thoroughly wet. Allow surface to remain undisturbed for at least 3 minutes before wiping. For use against HIV-1 see special instructions for cleaning.



## Disinfectant Aerosol and DEODORIZER

A BROAD-SPECTRUM DISINFECTANT FOR  
ENVIRONMENTAL SURFACES.

ELIMINATES ODORS BY ATTACKING  
THEM AT THEIR SOURCE.

DISINFECTANT • TUBERCULOCIDAL • \*VIRUCIDAL

DISINFECTANTE: SI USTED NO ENTIENDE LA ETIQUETA, BUSQUE A ALGUIEN PARA QUE SE LA EXPLIQUE A USTED EN DETALLE. (IF YOU DO NOT UNDERSTAND THIS LABEL, FIND SOMEONE TO EXPLAIN IT TO YOU IN DETAIL)

**ACTIVE INGREDIENTS:**

Orthophenylphenol:.....0.08%  
**INERT INGREDIENTS:..... 99.92%**  
**TOTAL:..... 100.00%**

### KEEP OUT OF REACH OF CHILDREN CAUTION

**CONTENTS UNDER PRESSURE**

SEE ADDITIONAL PRECAUTIONARY STATEMENTS ON SIDE PANEL

EPA REG. NO. 3862-188-82278

EPA EST. NO. 9143-TX-1

NET WT: 17 OZ (1 LB 1 OZ) (482 G)



Sold By: Bio-Tek

Division of ABC Compounding Co., Inc.

PO Box 16247 • Atlanta, GA 30321-0247 • Toll Free 800-795-9222

### STORAGE AND DISPOSAL

**GENERAL:** Do not contaminate other materials, including water, food or feed during transport, use, storage, or disposal. If damaged or leaking, soak up and wrap waste, then dispose of as below.

**Storage:** Keep in a cool, dry, locked area inaccessible to children and away from flames, sparks, sunlight, or other sources of heat or ignition. Containers heated above 130°F (54°C) may burst. **Disposal:** Refasten protective cap on container and discard in trash or offer for recycling, if available. Do not incinerate or puncture.

### PRECAUTIONARY STATEMENTS HAZARDS TO HUMANS AND DOMESTIC ANIMALS

**CAUTION.** Harmful if swallowed. Harmful if absorbed through skin. Causes moderate eye irritation. Avoid contact with skin, eyes or clothing. Wear safety gloves and glasses. Do not spray on or toward eyes or face. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. Remove and wash contaminated clothing before reuse. Do not contaminate food or feed, nor surfaces that they may contact.

FIRST AID:	
<b>If swallowed</b>	<ul style="list-style-type: none"> <li>• Call a poison control center or doctor for further treatment advice.</li> <li>• Have person sip a glass of water if able to swallow.</li> <li>• Do not induce vomiting unless told to do so by the poison control center or doctor.</li> <li>• Do not give anything by mouth to an unconscious person.</li> </ul>
<b>If on skin or clothing</b>	<ul style="list-style-type: none"> <li>• Take off contaminated clothing.</li> <li>• Rinse skin immediately with plenty of water for 15-20 minutes.</li> <li>• Call a poison control center or doctor for further treatment advice.</li> </ul>
<b>If in eyes</b>	<ul style="list-style-type: none"> <li>• Hold eye open and rinse slowly and gently with water for 15-20 minutes.</li> <li>• Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye.</li> <li>• Call a poison control center or doctor for further treatment advice.</li> </ul>
Have the product container or label with you when calling a poison control center or doctor, or going for treatment.	

### PHYSICAL AND CHEMICAL HAZARDS

Flammable. Contents under pressure. Keep away from heat, sparks and open flame. Do not puncture or incinerate container. Exposure to temperatures above 130°F may cause bursting. Do not use this product in or on electrical equipment due to the possibility of shock hazard.